

WHAT'S IN THE BOX

Portable Ultrasonic Nebulizer



Check that you have everything

The children's model and the adult model come with the same parts.

- ☐ **The nebulizer**
with its medication cup and the cap that covers the mist outlet
- ☐ **Adult mask**
the larger mask
- ☐ **Child mask**
the smaller mask
- ☐ **Mouthpiece**
for breathing in through the mouth
- ☐ **USB charging cable**
USB-A to USB-C
- ☐ **Instruction manual**



Children's model box



Adult model box



Side of the children's box: adult mask, child mask, and mouthpiece

Charge the nebulizer before its first use. The step-by-step guide, the video, and the original manuals are at wellness-supplies.com/how-to-use/ultrasonic.

Which model do you have?



Children's, blue



Children's, green



Children's, orange



Children's, white



Adult, orange